



medicalbees
System-driven Surgery
Operieren mit System

CF-MCC Carbon/PEEK

medical bees offers in addition to the well-known instrument sets made of stainless steel and titanium new instrument sets in PEEK/ Carbon for complete X-ray transparency and MRI-compatibility. CF stands for carbon fiber reinforced instrument material and thus makes the instrument resistant and safe.

medical bees bietet zusätzlich zu den bekannten Instrumentensets aus Edelstahl und Titan jetzt auch Instrumente aus PEEK/ Carbon für komplette Röntgentransparenz und MRT-Kompatibilität an. CF steht für kohlefaserverstärkte Kunststoffe und Carbonmaterialien, die das Instrument besonders robust und sicher machen.

- ✓ No artefacts in X-rays, CT-MRI-compatible
- ✓ Low weight combined with an excellent rigidity
- ✓ Excellent bio-compatibility and sterilization characteristics

- ✓ Komplette Röntgentransparenz und CT-/MRT-Kompatibilität
- ✓ Geringes Gewicht kombiniert mit einer hervorragenden Festigkeit und Dimensionsstabilität
- ✓ Exzellente Bio-Kompatibilität und Sterilisationseigenschaften

CARBON / PEEK RETRACTOR / RETRAKTOR

	Max. Spreizweite <i>max. opening</i>	Besonderheit <i>features</i>
MB1-10-997-15	85 mm 3,3"	knickbare Arme <i>arms with joints</i>

PEEK BLADES / PEEK VALVEN

	Länge x Breite <i>length x width</i>		Länge x Breite <i>length x width</i>
MB1-10-978-10	40 x 20 mm	MB1-10-978-55*	30 x 27 mm
MB1-10-978-15	50 x 20 mm	MB1-10-978-60	40 x 27 mm
MB1-10-978-20	60 x 20 mm	MB1-10-978-65	50 x 27 mm
MB1-10-978-25	70 x 20 mm	MB1-10-978-70	60 x 27 mm
		MB1-10-978-75	70 x 27 mm

*nicht im CF-MCC Set enthalten | *not included in CF-MCC Set!*

CF-MCC SET / TRAY

MB1-10-977-20	bestehend aus Tray MB1-20-968-80, Retractor MB1-10-997-15 und je 2 Valven <i>includes Tray MB1-10-968-80, Retractor MB1-10-997-15 and 2 blades each</i>
MB1-10-968-80	Siebkorb mit Deckel für CF-MCC Set <i>Tray with lid for CF-MCC Set</i> 243 x 253 x 50mm

