



# medicalbees<sup>🐝</sup> PRIME MALLETS



Neu:

- ✓ Erstmals Garantie auf Schweißnaht: wir geben 6 Jahre Garantie auf die neue Schweißnaht. Es handelt sich erstmals um eine unsichtbare, 4-fach verstärkte Schweißnaht.
- ✓ Bessere Ergonomie durch Hohlgriffdesign und damit verbunden eine Gewichtsreduzierung des Griffs um 20%.
- ✓ Bessere Reinigbarkeit durch Griffmulden statt Rändelgriff.
- ✓ Matt gebürstete Oberfläche.

New:

- ✓ For the first time warranty on weld: we give a 6 year warranty on the new weld - an invisible, 4-fold reinforced weld.
- ✓ Better ergonomics due to hollow grip design and associated weight reduction of the handle by 20%.
- ✓ Better cleaning through non-slip grip design instead of knurled handle.
- ✓ Brushed matt finish.

Art.Nr./Art.No.	Länge/Length	Breite/Width	Kopfgewicht/head weight
MB1-07-002-07	265mm	35mm	450g
MB1-07-002-08	265mm	50mm	900g

Stahlgriff mit gedrehten Griffmulden/ handle with special recessed grip.

Es stehen Ihnen stehen weitere 80 Hämmer zur Verfügung.  
Over 80 different mallet designs are available.