




LD-Set Caspar Micro Lumbar Discectomy Set



SST Stainless Steel

Produktinformationen Product information

LD-Set Caspar Micro Lumbar Discectomy Set

- Das LD-Set wurde speziell für die Aufnahme des Caspar-Instrumentariums zur Verwendung bei lumbaralen Diskektomien entwickelt - alle Instrumente in einem Set - zur Retraktion der posterioren Weichteile.
- Ein gefensterter Speculum für einen sicheren Sitz in der Muskulatur und ein Gegenretractor mit einem speziellen Snap-In Verschluss für eine sichere Verbindung der lateralen Valven stehen zur Verfügung.
- Zur Verkleinerung der Inzision und zur Begrenzung des Traumas für das Gewebe.


LD-Set Caspar Micro Lumbar Discectomy Set

- The LD-Set was specifically designed to accommodate the Caspar instrument offering for use during lumbar discectomy procedures – all instruments in one set – to retract posterior soft tissues
- A non-slip fenestrated specula for a secure fit in muscles and a counter retractor with a special snap closure for a secure connection of the lateral blade.
- To reduce the size of the incision and limit trauma to the tissue.

LD-SET
LD-SET

MB1-10-772-00

bestehend aus Produkten mit „*“ gekennzeichnet (siehe Rückseite)
consisting of products marked „*“ (see reverse page)



Product Special



CASPAR VALVE LATERAL
Caspar blade lateral

Artikelnummer <i>Item number</i>	Tiefe <i>depth</i>
MB1-10-774-10 *	40 mm
MB1-10-774-15 *	45mm
MB1-10-774-20 *	50 mm
MB1-10-774-25 *	55 mm
MB1-10-774-30 *	60 mm
MB1-10-774-35 *	65 mm
MB1-10-774-40 *	75 mm
MB1-10-774-45	80 mm
MB1-10-774-50	85 mm
MB1-10-774-55	90 mm
MB1-10-774-60	100 mm
MB1-10-774-65	110 mm
MB1-10-774-70	120 mm

schwarz, TiAIN beschichtet
black, TiAIN coated surface

CASPAR VALVE MEDIAL
Caspar blade medial

Artikelnummer <i>Item number</i>	Tiefe <i>depth</i>
MB1-10-773-05	40 mm
MB1-10-773-10	45mm
MB1-10-773-15	50 mm
MB1-10-773-20	55 mm
MB1-10-773-25	60 mm
MB1-10-773-30	65 mm
MB1-10-773-35	70 mm
MB1-10-773-40	75 mm


schwarz, TiAIN beschichtet
black, TiAIN coated surface

* included in Set MB1-10-772-00

made in GERMANY

Änderungen der technischen Angaben und Maße liegen in unserem Ermessen. Modelländerungen vorbehalten. Ohne unsere ausdrückliche schriftliche Genehmigung darf dieses Printmedium weder im Gesamten noch auszugsweise kopiert oder nachgedruckt werden.

Changes in technical specifications and dimensions are at our discretion. Subject to modification. Without our explicit written approval this print media may not be copied or reprinted as a whole or in extracts.



CASPAR

MB1-10-775-00 *

Gegensperrer, schwarz-silber (ohne Valven)
Counter retractor, black-silver (without blades)

CASPAR Wirbelsäulenspeculum
CASPAR spine specula black

Artikelnummer <i>Item number</i>	Blattlänge <i>blade length</i>
MB1-10-772-05 *	40 mm
MB1-10-772-10 *	45 mm
MB1-10-772-15 *	50 mm
MB1-10-772-20 *	55 mm
MB1-10-772-25 *	60 mm
MB1-10-772-30 *	65 mm
MB1-10-772-47	70 mm
MB1-10-772-35 *	75 mm
MB1-10-772-40	80 mm
MB1-10-772-45	85 mm
MB1-10-772-50	90 mm
MB1-10-772-55	100 mm
MB1-10-772-60	110 mm
MB1-10-772-65	120 mm

schwarz, TiAIN beschichtet
black, TiAIN coated surface



MB1-10-775-12 *

Entnahmepinzette für Valven
Blade extraction forceps

MB1-20-968-85 *

LD-Siebkorb mit Deckel, Maße: 537x250x100 mm (LxBxT), inkl. 2
Einsätze (ohne Instrumente)

*LD-Tray with lid, dimensions: 537x250x100 mm (LxWxD), incl. 2
inserts (without instruments)*

