

Product Special

KC-Retractor Set KC-Reaktor Set



Product information Produktinformationen

KC-Retractor Set

- For posterior lumbar and thoracic spine surgery.
- Standard Set with KC-Retractor and standard blades in different designs.
- 1 to 5 teeth design, various widths and depths of blades available.
- Optimum assembly, repositioning and stability.
- All instruments available in one tray, securely stored (dimensions: 540 x 253 x 50mm).

New:

- **KC-Lightweight Set with newly designed retractor** and connection of blades with additional safeguard.
- Retractor is 25% lighter compared to standard retractor due to new design.
- Blades feature an additional safeguard to allow fixation of blades.

made in **GERMANY**

Änderungen der technischen Angaben und Maße liegen in unserem Ermessen. Modelländerungen vorbehalten. Ohne unsere ausdrückliche schriftliche Genehmigung darf dieses Printmedium weder im Gesamten noch auszugsweise kopiert oder nachgedruckt werden.

Changes in technical specifications and dimensions are at our discretion. Subject to modification. Without our explicit written approval this print media may not be copied or reprinted as a whole or in extracts.

STANDARD KC-RETRACTOR SET

MB1-10-762-05	without tray, consisting of:
MB1-10-762-06	with tray and lid MB1-20-968-78, consisting of:

	type	width	length	quantity
MB1-10-762-10	Standard KC-Retractor			1
MB1-10-764-15	5 teeth	76 mm	64 mm	2
MB1-10-764-03	left	31 mm	38 mm	1
MB1-10-764-06	right	31 mm	38 mm	1
MB1-10-764-09	left	57 mm	44 mm	1
MB1-10-764-12	right	57 mm	44 mm	1

KC-Blades standard

LIGHTWEIGHT KC-RETRACTOR SET

MB1-10-762-13	without tray, consisting of:
MB1-10-762-14	with tray and lid MB1-20-968-78, consisting of:

	type	width	length	quantity
MB1-10-762-11	Lightweight KC-Retractor			1
MB1-10-764-65	5 teeth	76 mm	64 mm	2
MB1-10-764-54	left	31 mm	38 mm	1
MB1-10-764-56	right	31 mm	38 mm	1
MB1-10-764-59	left	57 mm	44 mm	1
MB1-10-764-62	right	57 mm	44 mm	1

KC-Blades with special safeguard

MB1-20-968-78	tray with lid, 540 x 253 x 50 mm
----------------------	----------------------------------

