

medical**bees**[®]
System-driven Surgery
Operieren mit System

Hämmer | Mallets





KIRK


190 mm – 7½"	Ø B	Kopfgewicht – Head weight
MB1-07-001-01	40 mm	630 g

Bronzekopf verchromt, massiv
Bronze head, chromium plated, massive



280 mm – 11"	Ø B	Kopfgewicht – Head weight
MB1-07-002-01	35 mm	450 g
MB1-07-002-02	50 mm	900 g

Stahlkopf, massiv
Stainless steel head, massive




185 mm – 7¼"	Ø B	Kopfgewicht – Head weight
MB1-07-002-03	40 mm	520 g

Stahlkopf, massiv
Stainless steel head, massive



265 mm – 10 1/2"	Ø B	Kopfgewicht – Head weight
MB1-07-002-05	35 mm	450 g
MB1-07-002-06	50 mm	900 g

Stahlkopf, massiv, mit Ferozellgriff
Stainless steel head, massive, with phenolic handle



BERGMANN

230 mm – 9"	Ø B	Kopfgewicht – Head weight
MB1-07-006-05	40 mm	260 g


Stahlkopf, massiv
Stainless steel head, massive



OMBREDANNE

265 mm – 10 1/2"	Ø B	Kopfgewicht – Head weight
MB1-07-006-07	40 mm	520 g

Stahlkopf, massiv
Stainless steel head, massive



265 mm – 10 1/2"	Ø A–B	Kopfgewicht – Head weight
MB1-07-010-01	42–30 mm	310 g

Stahlkopf, massiv
Stainless steel head, massive



265 mm – 10 1/2"	Ø B–A	Kopfgewicht – Head weight
MB1-07-010-05	42–30 mm	310 g

Stahlkopf, massiv, mit Ferozellgriff
Stainless steel head, massive, with phenolic handle



265 mm – 10 1/2"	Ø B	Kopfgewicht – Head weight
MB1-07-011-05	30 mm	290 g

Kopf mit Bleifüllung, rückschlagfrei
Lead filled head, reperussion free

MB1-07-011-20

Ø 30 (M 10) x 91 mm

Kunststoffbacke, alleine
Plastic disc, only



172 mm – 6¾"	Ø B	Kopfgewicht – Head weight
MB1-07-004-01	20 mm	100 g

Kupferkopf, massiv
Cooper head, massive



PARTSCH


180 mm – 7"	Ø B	Kopfgewicht – Head weight
MB1-07-013-10	22 mm	120 g

Stahlkopf, massiv
Stainless steel head, massive



172 mm – 6¾"	Ø B	Kopfgewicht – Head weight
MB1-07-004-05	20 mm	90 g

Stahlkopf, massiv
Stainless steel head, massive



MEAD

170 mm – 6¾"	Ø B	Kopfgewicht – Head weight
MB1-07-014-01	26 mm	180 g

Stahlkopf, auswechselbare Ferozellbacken
Stainless steel head, exchangeable discs of artificial wood

MB1-07-014-55

Ø 26 mm

Ferozellbacke, alleine
Ferozell faces, only



LUCAE

190 mm – 7½"	Ø B	Kopfgewicht – Head weight
MB1-07-012-01	19 mm	80 g
MB1-07-012-02	25 mm	140 g

Stahlkopf, massiv
Stainless steel head, massive



204 mm – 8"	Ø B	Kopfgewicht – Head weight
MB1-07-020-01	7 mm	80 g


Stahlkopf, massiv
Stainless steel head, massive



USA MODELL

183 mm – 7¼"	Ø B	Kopfgewicht – Head weight
MB1-07-030-01	30 mm	220 g

Stahlkopf, massiv
Stainless steel head, massive



USA MODELL

193 mm – 7¾"	Ø B	Kopfgewicht – Head weight
MB1-07-030-10	38 mm	900 g

Stahlkopf, massiv
Stainless steel head, massive



CLOWARD

191 mm – 7½"	Ø B	Kopfgewicht – Head weight
MB1-07-035-01	20 mm	120 g

Stahlkopf, massiv
Stainless steel head, massive



CLOWARD

229 mm – 9"	Ø B	Kopfgewicht – Head weight
MB1-07-035-05	25 mm	220 g

Stahlkopf, massiv
Stainless steel head, massive



COTTLE

180 mm – 7"	Ø B	Kopfgewicht – Head weight
MB1-07-040-02	30 mm	240 g


Stahlkopf, massiv, Aluminiumgriff, eloxiert
Stainless steel head, massive, Aluminium handle, anodized



COLLIN

204 mm – 8"	Ø B	Kopfgewicht – Head weight
MB1-07-022-01	30 mm	150 g

Stahlkopf, massiv
Stainless steel head, massive




227 mm – 9"	Ø B	Kopfgewicht – Head weight
MB1-07-008-01	32 mm	140 g

Eloxiertes Aluminiumgriff und Aluminiumkopf
Anodized aluminium handle and head

MB1-07-008-02	
Ø 32 (M 10 mm)	20 g

Kunststoffbacken, alleine
Plastic discs, only



180 mm – 7"	Ø B	Kopfgewicht – Head weight
MB1-07-025-01	25 mm	130 g

Stahlkopf, massiv, 2 austauschbare Kunststoffbacken
Aluminiumgriff, eloxiert

Stainless steel head, massive, 2 exchangeable plastic discs
Aluminium handle, anodized

MB1-07-025-50	
Ø 25 (M 20 x 1 mm)	10 g

Kunststoffbacken, alleine
Plastic discs, only



185 mm – 7 1/4"	Ø B	Kopfgewicht – Head weight
MB1-07-024-01	25 mm	130 g

Stahlkopf, massiv
Stainless steel head, massive

MB1-07-024-50	
Ø 25 (M 20 x 1 mm)	10 g

Kunststoffbacke, alleine
Plastic disc, only



180 mm – 7"	Ø B	Kopfgewicht – Head weight
MB1-07-025-05	25 mm	140 g

Stahlkopf, massiv, 1 austauschbare Kunststoffbacke
Aluminiumgriff, eloxiert

Stainless steel head, massive, 1 exchangeable plastic disc
Aluminium handle, anodized

MB1-07-025-50	
Ø 25 (M 20 x 1 mm)	10 g

Kunststoffbacke, alleine
Plastic disc, only



185 mm – 7 1/4"	Ø B	Kopfgewicht – Head weight
MB1-07-024-05	38 mm	470 g

Stahlkopf, massiv
Stainless steel head, massive

MB1-07-024-55	
Ø 38 (M 24 x 1,5 mm)	30 g

Kunststoffbacke, alleine
Plastic disc, only



		Ø B	Kopfgewicht – Head weight
MB1-07-024-10	195 mm – 7¾"	25 mm	140 g
MB1-07-024-15	200 mm – 8"	32 mm	250 g

Stahlkopf, massiv, Griff hohl
Stainless steel head, massive, hollow handle

	Ø B	Kopfgewicht – Head weight
MB1-07-024-40	25 (M 18 x 1 mm)	10 g
MB1-07-024-45	32 (M 18 x 1 mm)	20 g

Kunststoffbacke, alleine
Plastic disc, only



195 mm – 7¾"	Ø B	Kopfgewicht – Head weight
MB1-07-018-10	25 mm	110 g

Stahlkopf, massiv, Griff hohl
Stainless steel head, massive, hollow handle

	Ø B	Kopfgewicht – Head weight
MB1-07-018-50	25 mm (M 20 x 1.5 mm)	10 g

Kunststoffbacke, alleine
Plastic disc, only