



Rippenspreizer | Rib Spreaders





BAILEY

MB1-10-700-06

180 mm – 7"

A = 19 mm

BAILEY GIBBON

MB1-10-700-11

180 mm – 7"

A = 29 mm

BAILEY PEDIATRIC BABY

MB1-10-701-10

140 mm – 5 ½ "

A = 18 mm



LEMMON STERNAL

MB1-10-701-18

145 mm – 5 ¾ "

Spreizer für Kinder
Retractor for children


MB1-10-701-20

180 mm – 7"

Groß
Large





Rippenspreizer | Rib Spreaders



TUFFIER

	S	L	B	T
MB1-10-704-05	178 mm	108 mm	45 mm	48 mm


Gekerbte Valven
Roughened blades



TUFFIER

	S	L	B	T
MB1-10-704-10	178 mm	107 mm	45 mm	48 mm


Geriefte Valven
Grooved blades



BADGLEY



	S	L	B	T
MB1-10-704-15	165 mm	108 mm	45 mm	50 mm

Valve mit vier Zähnen
Blade with four prongs




RIENHOFF

	S	L	B	T
MB1-10-705-05	140 mm	116 mm	33 mm	30 mm




RIENHOFF

	S	L	B	T
MB1-10-705-10	140 mm	146 mm	45 mm	47 mm



HAIGHT

	S	L	B	T
MB1-10-710-05	126 mm	64 mm	31 mm	32 mm
MB1-10-710-10	126 mm	90 mm	31 mm	32 mm

**BURFORD**

	S	L	B	T
MB1-10-731-12	208 mm	192 mm		

Inklusive je 1 Paar Seitenvalven
Including 1 pair each of blades


MB1-10-731-41			65 mm	77 mm
MB1-10-731-51			65 mm	46 mm

BURFORD

	S	L	B	T
MB1-10-734-05	200 mm	205 mm		

Inklusive je 1 Paar gelochter Seitenvalven, schräg, mit 2 Spitzen
Including 1 pair each of perforated blades, inclined, with 2 spikes

MB1-10-734-50			64 mm	65 mm
MB1-10-734-60			64 mm	80 mm

**BURFORD-FINOCHIETTO CHILD**

	S	L	B	T
MB1-10-731-70	157 mm	108 mm		

Inklusive je 1 Paar Seitenvalven
Including 1 pair each of blades

MB1-10-731-80			44 mm	47 mm
MB1-10-731-85			45 mm	31 mm
MB1-10-731-90			33 mm	30 mm
MB1-10-731-95			23 mm	22 mm





MB1-10-731-90



MB1-10-731-85



MB1-10-731-95



FINOCHIETTO BABY

	S	L	B	T
MB1-10-715-01	55 mm	51 mm	15 mm	13 mm
MB1-10-715-03	75 mm	51 mm	15 mm	13 mm
MB1-10-715-05	80 mm	60 mm	21 mm	18 mm

Gerade Seitenarme
Straight arms

	S	L	B	T
MB1-10-715-07	80 mm	60 mm	21 mm	18 mm

Gebogene Seitenarme
Curved arms

FINOCHIETTO CHILD



	S	L	B	T
MB1-10-717-50	150 mm	111 mm	32 mm	28 mm

Gebogene Seitenarme, mit rundem Fenster
Curved arms, with circular fenestration




FINOCHIETTO BABY

	S	L	B	T
MB1-10-717-05	80 mm	68 mm	21 mm	19 mm

**FINOCHIETTO**



	S	L	B	T
MB1-10-718-10	200 mm	186 mm	66 mm	51 mm
MB1-10-718-15	265 mm	186 mm	66 mm	51 mm
MB1-10-718-20	310 mm	186 mm	66 mm	51 mm

Gerade Seitenarme
Straight arms

**FINOCHIETTO**


	S	L	B	T
MB1-10-718-06	148 mm	130 mm	45 mm	40 mm
MB1-10-718-07	148 mm	136 mm	45 mm	40 mm
MB1-10-718-50	200 mm	162 mm	65 mm	50 mm
MB1-10-718-60	200 mm	130 mm	65 mm	50 mm
MB1-10-718-25	200 mm	180 mm	65 mm	50 mm
MB1-10-718-55	220 mm	180 mm	65 mm	70 mm
MB1-10-718-57	250 mm	180 mm	65 mm	70 mm
MB1-10-718-30	265 mm	180 mm	65 mm	50 mm
MB1-10-718-35	310 mm	180 mm	65 mm	50 mm

Gebogene Seitenarme
Curved arms



MORSE CHILD

	S	L	B	T
MB1-10-720-01	154 mm	125 mm	22 mm	15 mm



MORSE


	S	L	B	T
MB1-10-720-05	212 mm	148 mm	29 mm	23 mm

FAVALORO-MORSE

	S	L	B	T
MB1-10-720-10	212 mm	148 mm	29 mm	38 mm
MB1-10-720-15	212 mm	205 mm	29 mm	38 mm

PACIFICO-MORSE

	S	L	B	T
MB1-10-720-20	170 mm	148 mm	29 mm 29 mm	23 mm 34 mm



DE BAKEY ORIGINAL

	S	L	B	T
MB1-10-722-00	192 mm	153 mm		

Inklusive je 1 Paar Seitenvalven
Including 1 pair each of blades

MB1-10-722-10			30 mm	30 mm
MB1-10-722-15			40 mm	40 mm
MB1-10-722-20			50 mm	50 mm




DE BAKEY ORIGINAL

	S	L	B	T
MB1-10-724-00	240 mm	170 mm		

Inklusive je 1 Paar Seitenvalven
Including 1 pair each of blades


MB1-10-724-10			60 mm	60 mm
MB1-10-724-15			80 mm	50 mm
MB1-10-724-20			100 mm	40 mm



ANKENY small

	S	L	B	T
MB1-10-706-07	139 mm	114 mm	19 mm	23 mm


Mit 4 fest montierten Wechselvalven, natur eloxiertes Aluminium
With 4 fixed swivel blades, clear anodized aluminium



ANKENY large

	S	L	B	T
MB1-10-706-12	195 mm	289 mm	35 mm	35 mm


Mit 6 drehbaren Valven
With 6 swivel blades



HAIGHT-FINOCHIETTO

	S	L	B	T
MB1-10-712-05	98 mm	76 mm	31 mm	31 mm

Natur eloxiertes Aluminium
Clear anodized aluminium


**BURFORD-FINOCHIETTO**

	S	L	B	T
MB1-10-732-05	135 mm	113 mm		

Inklusive je 1 Paar Seitenvalven
Including 1 pair each of blades

MB1-10-732-30			42 mm	34 mm
MB1-10-732-40			42 mm	48 mm

natur eloxiertes Aluminium
clear anodized aluminium


**BURFORD-FINOCHIETTO**

	S	L	B	T
MB1-10-732-10	200 mm	197 mm		
MB1-10-732-15	288 mm	197 mm		
MB1-10-732-20	330 mm	197 mm		

Inklusive je 1 Paar Seitenvalven
Including 1 pair each of blades


MB1-10-732-50			62 mm	47 mm
MB1-10-732-60			62 mm	65 mm

natur eloxiertes Aluminium
clear anodized aluminium



CASTANEDA

	S	L	B	T
MB1-10-736-05	74 mm	60 mm	34 mm	11 mm
MB1-10-736-10	80 mm	80 mm	45 mm	13 mm
MB1-10-736-15	110 mm	126 mm	60 mm	17.5 mm




COOLEY CHILD

	S	L	B	T
MB1-10-738-05	74 mm	58 mm	17 mm	11 mm
MB1-10-738-10	101 mm	80 mm	30 mm	17 mm
MB1-10-738-15	74 mm	58 mm	17 mm	13 mm



COOLEY CHILD



	S	L	B	T
MB1-10-738-20	97 mm	77 mm	20.5 mm	17 mm



BAILEY FINOCHIETTO

	S	L	B	T
MB1-10-738-50	74 mm	58 mm	14 mm	12 mm

natur eloxiertes Aluminium
clear anodized aluminium



COOLEY STERNUM

	S	L	B	T
MB1-10-708-06	137 mm	112 mm	46 mm	25 mm
MB1-10-708-10	170 mm	140 mm	52 mm	34 mm
MB1-10-708-15	194 mm	158 mm	102 mm	29 mm





DE BAKEY klein

	S	L	B	T
MB1-10-721-05	152 mm	144 mm		

Inklusive je 1 Paar Seitenventilen
Including 1 pair each of blades

MB1-10-721-30			30 mm	30 mm
MB1-10-721-40			40 mm	40 mm
MB1-10-721-50			50 mm	50 mm



DE BAKEY groß


	S	L	B	T
MB1-10-723-05	241 mm	189 mm		

Inklusive je 1 Paar Seitenvalven
Including 1 pair each of blades

MB1-10-723-60			60 mm	60 mm
MB1-10-723-50			80 mm	50 mm
MB1-10-723-40			100 mm	40 mm

optional
optional

MB1-10-723-80			60 mm	80 mm
MB1-10-723-90			66 mm	50 mm



MB1-10-723-50




MB1-10-723-40



MB1-10-723-80



MB1-10-723-90





HARKEN Rippen- und Scapulaspreizer / Rib and scapula spreader

	S	L	B	T
MB1-10-728-05	280 mm	210 mm	65 mm	49 mm
Inklusive 1 Paar abnehmbarer Valven Including 1 pair each of detachable blades				
MB1-10-728-15			50 mm	47 mm



DUBOST Thoraxsperrer

	S	L	B	T
MB1-10-725-30	160 mm	140 mm		
Inklusive je 1 Paar auswechselbaren Valven Including 1 pair each of blades				
MB1-10-725-40			100 mm	30 mm
MB1-10-725-50			100 mm	40 mm
MB1-10-725-60			100 mm	50 mm



MB1-10-727-10



MB1-10-727-50

DUBOST Thoraxsperrer



	S	L	B	T
MB1-10-725-80	160 mm	140 mm		
Inklusive je 1 Paar Seitenvalven Including 1 pair each of blades				
MB1-10-725-40			100 mm	30 mm
MB1-10-725-50			100 mm	40 mm
MB1-10-725-60			100 mm	50 mm

Modell inkl. Mammariavalve (MB1-10-727-10) und Aufnahme (MB1-10-727-50)
Model incl. mammary blade (MB1-10-727-10) and fixation blade (MB1-10-727-50)







Rippenspreizer | Rib Spreaders



Mammaria-Valve
Mammary blade

IMA B-RETRACTOR

	S	L	B	T
MB1-10-727-05	134 mm	195 mm	150 mm	50mm



IMA RETRACTOR

MB1-10-729-06

komplett, Spreizweite: 230 mm
inklusive

- 1 Rahmen mit Armen (MB1-10-729-04)
- 1 Doppelaufsteckvalve (MB1-10-729-09)
- 1 Doppel-Aufsteckvalve mit Nute zur Befestigung des T-Gestells (MB1-10-729-07)
- 1 T-Stück (MB1-10-729-15)
- 2 Universalklemmen mit Verzahnung (MB1-10-729-30)
- 2 COOLEY Drahtvalven klein, 25 mm (MB1-10-729-35)
- 1 COOLEY Drahtvalve medium, 40 mm (MB1-10-729-40)
- 1 COOLEY Drahtvalve groß, 50 mm (MB1-10-729-45)
- 2 biegsame RIBBON Valven (MB1-10-729-50)

complete, spreading width: 230 mm

including:

- 1 Frame with arms (MB1-10-729-04)
- 1 Double fixation blade (MB1-10-729-09)
- 1 Double fixation blade with notch for fixation to T-rack (MB1-10-729-07)
- 1 T-clamp (MB1-10-729-15)
- 2 Universal clamps (MB1-10-729-30)
- 2 COOLEY small blade retractors, 25 mm (MB1-10-729-35)
- 1 COOLEY medium blade retractors, 40 mm (MB1-10-729-40)
- 1 COOLEY large blade retractors, 50 mm (MB1-10-729-45)
- 2 RIBBON malleable blades (MB1-10-729-50)



